



**METALIC GOLD RINGING GEL
ETB21125**

INGREDIENTS	INCI NAME	
Water	Water	55.0%
Glycerin	Glycerin	10.0%
Drakeol 10 light	Mineral oil	12.0%
Prespearl Deep Gold	Mica & Titanium Dioxide & Iron Oxide	1.0%
Hetoxol OL-2	Oleth-2	4.4%
Hetoxol I-20-20	Octyldodeceth-20	<u>17.6%</u>
		100.0%

PROCEDURE

- 1) Combine the water and Glycerin and heat to 80C.
- 2) In a separate vessel combine the Mineral Oil, Hetoxol OL-2 and Hetoxol I-20-20 and heat to the same temperature.
- 3) With adequate mixing add step 2 to step 1 and maintain at 80C for 5 minutes.
- 4) Reduce mixing speed and cool to 40C, add Prespearl Deep Gold.
- 5) Cool to room temperature.

Global Seven, Inc. believe the information contained in this literature to be reliable, but it is only representative information. No warranty is given or implied as to accuracy, suitability for particular applications or the results to be obtained. This is not a commercialized product and has not been the subject of extensive testing for efficacy, stability, performance or safety. No recommendation should be construed as inducement to infringe patents.