



**MOISTURIZING ANTI-BACTERIAL  
LOTION with TRICLOSAN  
ETB567**

<b>INGREDIENTS</b>	<b>INCI NAME</b>	
Water	Water	24.65%
Ultrez 21	Acrylates/C10-30 Alkyl Acrylate Crosspolymer	0.80%
SD Alcohol 40-B	Alcohol Denat.(Ethyl Alcohol)	63.00%
Triclosan	Triclosan	1.00%
AMP-95	Aminomethyl Propanol	0.55%
Hest IS-3-O	Isosteareth-3 Ethylhexanoate	<u>10.00%</u>
		100.00%

**PROCEDURE**

- 1) Add the Ultrez-21 to water and once hydrated begin mixing slowly.
- 2) With slow mixing add the Ethanol and mix until uniform.
- 3) Add the Triclosan and mix until uniform.
- 4) Add the AMP-95 and mix until uniform.
- 5) Add the Hest IS-3-O and mix slowly until uniform.

Global Seven, Inc. believes the information contained in this literature to be reliable, but it is only representative information. No warranty is given or implied as to accuracy, suitability for particular applications or the results to be obtained. This is not a commercialized product and has not been the subject of extensive testing for efficacy, stability, performance or safety. No recommendation should be construed as inducement to infringe patents.